

# Swing Away SYSTEM




P17W656R2



## Swing Away

Swing Away Features.....	4
Facial Items Adjustment.....	5
Swing Away Accessories and Options .....	7
Functionality Meets Convenience.....	8
Special Situations .....	9



# Swing Away SYSTEM

**Helps you place pads,  
switches, sensor  
switches, sip & puff,  
alternative driving  
methods, or lateral support  
where it is needed most**

SAR2B-B shown with the All Positioning (AP) head pad AP350M. Installed in the left side rod a 910M Spot Facial Pad and on the right side rod a red Egg Switch ESRED.



## **SAR2A Series**

2 points standard profile Swing Away is compatible with the Comfort Plus Contour (CP).



## **SAR2B Series**

The 2 points low profile Swing Away is compatible with the All Positioning (AP), Combo (CB), i2i, and UniLink Ultra head support.



## **SAR2C Series**


The 2 points standard flat profile Swing Away is compatible with the Comfort Plus Series Adjustable (ACP) and the Tri-Comfort with Fixed Angled Laterals (TCP).



## **SAR4A Series**

4 points profile Swing Away is compatible with the Comfort Plus Contour (CP) head support.

# Swing Away Features



**SAR2A Plate**  
Standard Profile

**SAR2B Plate**  
Low Profile

**SAR4A Plate**  
4pt Profile



Two steps  
30° increments  
of height adjustment



# Facial Items Adjustment

SAR2B-B is included in the i-Drive IDHBT300-U as shown. The 17mm ball link is included on the head support hardware and is not included with the i-Drive head array.




Our cord wrap keeps switch wires tight against the rod, reducing risk of these being caught and damaged. The wrap included with the i-Drive is made out of Neoprene.



**120°**

APPROXIMATE  
FREE ANGLE FROM  
ALL DIRECTIONS

SAR2B-B is included in the i-Drive IDHBT300-U as shown. The 17mm ball link is included on the head support hardware and is not included with the i-Drive head array.



**92°**

Easy to Use  
Companion  
Assist Release  
Mechanism

# Swing Away Accessories and Options

Stealth Products offers a series of mounting options for facial pads, sensors, joysticks and sip & puff mounting solutions. Here are some examples.



**SP265M**

Sip & Puff extension for the SAR system

**IDH100-48M**

Facial Pad with Embedded Proximity Sensor Switch for i-Drive®

**F8675M**

Figure 8 Egg Switch Mount allows for a second rod to be added to one side.

Red Egg Switch

**ESRED**

## Facial Pads



**930M 920M 900M 910M**



For lateral or anterior head support

## Rod Options

Rods are available from 6in (15.24cm) to 13in (33.02cm) in length.

SU262 **6in** (15.24cm)

SU260 **9in** (22.86cm)

SU264 **13in** (33.02cm)

Extended Mounting Rods are available in 2 lengths. Mostly used for joysticks.

SU265 **18in** (45.72cm)

SU266 **24in** (60.96cm)




Also ask for our wide variety of switches, sensors, joystick and i-Drive®

# Functionality Meets Convenience

The Swing Away system is used with the i-Drive Pro Sip & Puff Head Array (IDHBT310-U) to mount both the rod and facial pad as well as the sip & puff flex tubing.

The swing away system provides functionality without compromising the user's field of view.



The new Swing Away system works with all Stealth Products head solutions. Compatible with 17mm ball head pad configurations.



# Special Situations

The Swing Away system allows every element to be precisely positioned to address even the most challenging situations.



Stealth Products has several mounting plates for the rods, including the Egg, Microlight (on flex tube), and other 3<sup>rd</sup> party switches, like AbleNet or MoVis. **Fiber Optics** can also be placed in the flex tube.

**Each side adjusts  
independently**

SAR2B-B shown with  
the i2i Torso Positioning  
System I2I555M.  
Installed in the left side  
rod a 910M Facial Pad  
and on the right side  
rod a red Egg Switch  
ESRED. Mounted on  
UniLink Flip-Down Head  
Solutions Hardware

• **Facial Elements  
can easily be  
positioned or rotated  
to the user's head  
orientation**





<https://stealthproducts.com>

+1 (800) 965-9229

+1 (512) 715-9995

[info@stealthproducts.com](mailto:info@stealthproducts.com)



**104 John Kelly Drive, Building A, Burnet, TX 78654**